



# BioBasic AX HPLC Columns

Optimized for the separation of proteins, peptides, other anionic species and polar molecules



- ▶ Weak anion exchange phase for multiple charged species
- ▶ 300Å pore size for enhanced protein and peptide separations
- ▶ Suitable for HILIC retention and separation of highly polar molecules
- ▶ Superb stability under demanding pH conditions



BioBasic AX HPLC Columns					
Length	1.0mm I.D.	2.1mm I.D.	3.0mm I.D.	4.0mm I.D.	4.6mm I.D.
<b>Particle Size 5μm</b>					
50mm	73105-051030	73105-052130	73105-053030	73105-054030	73105-054630
100mm	73105-101030	73105-102130	73105-103030	73105-104030	73105-104630
150mm	73105-151030	73105-152130	73105-153030	73105-154030	73105-154630
250mm	73105-251030	73105-252130	73105-253030	73105-254030	73105-254630

Other column dimensions are available. Please call Customer Service for more information.


BioBasic AX Drop-in Guard Cartridges						
Particle Size	Length	4.6mm/4.0mm ID	3.0mm ID	2.1mm ID	1.0mm ID	Quantity
5μm	10mm	73105-014001	73105-013001	73105-012101	73105-011001	4 Pack
		850-00	852-00	852-00	851-00	1 Each




BioBasic AX KAPPA Capillary HPLC Columns					
Length (mm)	75μm I.D.	100μm I.D.	180μm I.D.	320μm I.D.	500μm I.D.
<b>Particle Size 5μm</b>					
100mm	73105-100065	73105-100165	73105-100265	73105-100365	73105-100565
150mm	73105-150065	73105-150165	73105-150265	73105-150365	73105-150565

## BioBasic SCX Columns

For the separation of proteins, peptides, and other cationic species



- ▶ Strong cation exchange phase based on sulfonic acid chemistry
- ▶ Separation and retention of basic and other cationic species
- ▶ 300Å pore size for enhanced protein and peptide separations
- ▶ Outstanding stability under demanding pH conditions



BioBasic SCX Cation Exchange HPLC Columns					
Length	1.0mm I.D.	2.1mm I.D.	3.0mm I.D.	4.0mm I.D.	4.6mm I.D.
<b>Particle Size 5µm</b>					
50mm	73205-051030	73205-052130	73205-053030	73205-054030	73205-054630
100mm	73205-101030	73205-102130	73205-103030	73205-104030	73205-104630
150mm	73205-151030	73205-152130	73205-153030	73205-154030	73205-154630
250mm	73205-251030	73205-252130	73205-253030	73205-254030	73205-254630

Other column dimensions are available. Please call Customer Service for more information.

BioBasic SCX Guard Cartridges						
Particle Size	Length	4.0mm ID	3.0mm ID	2.1mm ID	1.0mm ID	Quantity
5µm	10mm	73205-014001	73205-013001	73205-012101	73205-011001	4 Pack
	UNIGUARD Drop-in Guard Cartridge Holder	850-00	852-00	852-00	851-00	1 Each

BioBasic SCX KAPPA HPLC Columns					
Length	75µm I.D.	100µm I.D.	180µm I.D.	320µm I.D.	500µm I.D.
<b>Particle Size 5µm</b>					
50mm	73205-050065	73205-050165	73205-050265	73205-050365	73205-050565
100mm	73205-100065	73205-100165	73205-100265	73205-100365	73205-100565
150mm	73205-150065	73205-150165	73205-150265	73205-150365	73205-150565

BioBasic SCX KAPPA Guard Columns			
Length	180µm I.D.	320µm I.D.	500µm I.D.
<b>Particle Size 5µm</b>			
30mm	73205-030215	73205-030315	73205-030515

BioBasic SCX Nanobore HPLC Columns							
Particle Size	Length	75µm ID	75µm ID Multipack	Quantity	150µm ID	150µm ID Multipack	Quantity
<b>IntegraFrit</b>							
5µm	50mm	73205-057563	73205-057564	3 Pack	73205-051563	73205-051564	3 Pack
5µm	100mm	73205-107563	73205-107564	3 Pack	73205-101563	73205-101564	3 Pack
<b>PicoFrit, 15µm ID Tip</b>							
5µm	50mm	73205-057581	73205-057582	3 Pack	--	--	--
5µm	100mm	73205-107581	73205-107582	3 Pack	--	--	--

Unless otherwise specified, IntegraFrit and PicoFrit are sold in single-column units.